# Package 'affyQCReport'

April 11, 2018		
Version 1.56.0		
Title QC Report Generation for affyBatch objects		
Author Craig Parman <craig.parman@bifx.org>, Conrad Halling <conrad.halling@bifx.org>, Robert Gentleman</conrad.halling@bifx.org></craig.parman@bifx.org>		
Maintainer Craig Parman <craig.parman@bifx.org></craig.parman@bifx.org>		
<b>Depends</b> Biobase (>= 1.13.16), affy, lattice		
<b>Imports</b> affy, affyPLM, Biobase, genefilter, graphics, grDevices, lattice, RColorBrewer, simpleaffy, stats, utils, xtable		
Suggests tkWidgets (>= 1.5.23), affydata (>= 1.4.1)		
<b>Description</b> This package creates a QC report for an AffyBatch object.  The report is intended to allow the user to quickly assess the quality of a set of arrays in an AffyBatch object.		
License LGPL (>= 2)		
biocViews Microarray,OneChannel,QualityControl		
Collate arrayReport.R correlation.R border1.R border2.R titlepage.R signaldist.R affyQAReport.R		
NeedsCompilation no		
R topics documented:		
affyQAReport		
borderQC1		
borderQC2		
correlationPlot		
QCReport		
signalDist		
titlePage		
Index		

2 affyQAReport

affyQAReport	Generate a QA Report for Affymetrix arrays	

#### **Description**

These functions provide a method for creating, viewing and removing QA reports on an experiment run on Affymetrix CEL files.

#### Usage

#### **Arguments**

affyB	An instance of the AffyBatch class.
output	Type of output wanted, currently only pdf is supported.
outdir	The directory (or folder) where output will be directed.

overwrite Whether an existing report, with the same name, should be overwritten.

repName A character string indicating the name of the report.

x A character string containing the location of the directory holding the report.

#### **Details**

A wide variety of tools from different packages are used, in conjunction with Sweave tools to produce an integrated document with various statistics that should help determine whether there are data quality problems with the data.

This package only works for Affymetrix data. It requires a working LaTeX implementation and so is not likely to work on Windows platforms, and in some cases on OS X unless the user has installed these tools.

Sample names that have an underscore in them have it replaced by a dot. Underscore causes a number of problems with Latex, if it is not escaped.

#### Value

A list with the following components is returned.

qcStats The QC statistics computed by the **simpleaffy** package. affyPLM The QC statistics computed by the **affyPLM** package.

MADS The computed MAD statistic on the 'M' values from an M-A decomposition.

loc The location (directory/folder) containing the report.

name The name of the report.

## Author(s)

R. Gentleman

borderQC1 3

## **Examples**

```
## Not run:
    library("affydata")
    data("Dilution")
    affyQAReport(Dilution)
## End(Not run)
```

borderQC1

Distribution of intensities of the border elements

## Description

Creates the forth page of the QC report created by QCReport for an object of class AffyBatch The page contains a boxplots showing the intensities from the border element of arrays in an AffyBatch object.

## Usage

```
borderQC1(object)
```

## Arguments

object

An object of class AffyBatch.

#### **Details**

This creates the forth page of a QCReport report. Boxplots of the positive and negative border elements of arrays in a AffyBatch object are shown.

More details are available in the vignette.

## Value

**TRUE** 

## **Examples**

```
library(affydata)
data(Dilution)
```

borderQC1(Dilution)

4 correlationPlot

borderQC2

Center of intensity QC plots

## Description

Creates the fifth page of the QC report created by QCReport for an object of class AffyBatch The page contains plots showing the center of intensity of the border element of arrays in an AffyBatch object.

## Usage

```
borderQC2(object)
```

#### **Arguments**

object

An object of class AffyBatch.

#### **Details**

This creates the fifth page of a QCReport report. Plots of the center of intensity of the positive and negative border elements of arrays in a AffyBatch object are shown. This plot is useful for detecting spatial biases in intensities on an array.

More details are available in the vignette.

#### Value

**TRUE** 

### **Examples**

```
library(affydata)
data(Dilution)
borderQC2(Dilution)
```

correlationPlot

Array-array intensity correlation plot

#### **Description**

Creates the sixth page of the QC report created by QCReport for an object of class AffyBatch The page displays array-array intensity correlation coefficients.

#### Usage

```
correlationPlot(object)
```

QCReport 5

#### **Arguments**

object An object of class AffyBatch.

#### **Details**

This creates the sixth page of a QCReport report. The page displays array-array intensity correlation coefficients for all pairs of arrays in an AffyBatch object. Arrays are ordered according to phenotypic data if available. This plot is useful for detecting outlier arrays.

More details are available in the vignette.

#### Value

**TRUE** 

#### **Examples**

```
library(affydata)

data(Dilution)

correlationPlot(Dilution)
```

QCReport

QC Report for AffyBatch objects

## Description

Creates a QC report for an object of class AffyBatch.

#### Usage

```
QCReport(object,file="AffyQCReport.pdf",...)
```

#### **Arguments**

 $object \hspace{1cm} An \hspace{1cm} object \hspace{1cm} of \hspace{1cm} class \hspace{1cm} \textbf{AffyBatch}.$ 

file A valid filename for the output PDF file.
... further arguments for plot function.

#### **Details**

This creates a 6 page QC report in PDF file format.

- The default output file is AffyQCReport.pdf created in the working directory.
- The AffyBatch object is not required to contain phenotypic data but that data is used if available.
- The CDF environment must be available for the method to work fully.

More details are available in the vignette.

6 signalDist

#### Value

**TRUE** 

## Examples

```
## Not run:
    library(affydata)
    data(Dilution)

QCReport(Dilution,file="ExampleQC.pdf") ##A QC report for the Example data will be generated
## End(Not run)
```

signalDist

Intensity distribution plots for a AffyBatch object

#### **Description**

Creates the secondpage of the QC report created by QCReport for an object of class AffyBatch The page contains a boxplot and intensity distribution plots.

#### Usage

```
signalDist(object)
```

#### **Arguments**

object

An object of class AffyBatch.

## **Details**

This creates the second page of a QCReport report. The first graph is created with the AffyBatch method boxplot. The second graph in created with the AffyBatch method hist.

More details are available in the vignette.

#### Value

**TRUE** 

## **Examples**

```
library(affydata)
data(Dilution)
signalDist(Dilution)
```

titlePage 7

titlePage

QC report title page with array names

## Description

Creates the titlepage of the QC report created by QCReport for an object of class AffyBatch.

## Usage

```
titlePage(object)
```

## Arguments

object

An object of class AffyBatch.

#### **Details**

This creates the title page of a  ${\tt QCReport}$  report.

More details are available in the vignette.

## Value

TRUE

## **Examples**

```
library(affydata)
data(Dilution)
```

titlePage(Dilution)

## **Index**

```
*Topic hplot
    borderQC1, 3
    borderQC2, 4
    correlationPlot, 4
    QCReport, 5
    \verb|signalDist|, 6
    titlePage, 7
*Topic manip
    {\it affyQAReport, 2}\\
AffyBatch, 3-7
{\it affyQAReport, 2}\\
borderQC1, 3
borderQC2, 4
{\tt correlationPlot}, 4
openQAReport (affyQAReport), 2
plot, 5
QCReport, 3–5, 5, 6, 7
QCReport-correlationPlot
         (correlationPlot), 4
QCReport-signalDist(signalDist), 6
QCReport-titlePage (titlePage), 7
rmQAReport (affyQAReport), 2
signalDist, 6
titlePage, 7
```