MAQCsubsetILM: MAQC reference subset for the Illumina platform

Laurent Gatto

November 3, 2022

1

 $\mathbf{2}$

Contents

1 The MAQC reference datasets

2 Loading the reference data

1 The MAQC reference datasets

The MAQC (MicroArray Quality Control) project¹ provides a set of reference datasets for a set of 10 platforms (see *Summary of the MAQC Data Sets*² for more details). This package provides a subset of the Illumina MAQC dataset³.

Regarding the Illumina platform (ILM prefix), a total of 59 Human-6 BeadChip 48K v1.0 have been generated. Four different reference RNAs have been used: (A) 100% of Stratagene's Universal Human Reference RNA, (B) 100% of Ambion's Human Brain Reference RNA, (C) 75% of A and 25% of B and (D) 25% of A and 75% of B. Each reference has been repeated 5^4 times (noted _A1_ to _A5_)⁵ on three different test sites (noted _1_ to _3_). As an example, the .CEL result file for the first replicate of test site 2, for the reference ARN C is named ILM_2_C1.CEL.

¹http://www.fda.gov/nctr/science/centers/toxicoinformatics/maqc

²http://edkb.fda.gov/MAQC/MainStudy/upload/Summary_MAQC_DataSets.pdf

³Packages for the datasets of other platforms will follow and will all be named MAQCsubsetXXX where XXX is the three-letter code used by the MAQC consortium.

⁴except for site 1, reference C, where 4 replicates are available

⁵the replicates for site 2, reference D are labelled $_D1_, _D2_, _D4_, _D6_$ and $_D7_$

These datasets are freely available and allow, for example, researchers to compare the reproducibility of their own Human-6 BeadChip 48K v1.0 data with the MAQC data. *MAQCsubsetILM* offers 3 randomly chosen BeadChips for each reference RNA, one for each test site. Each reference RNA subset is accessible as an R data object, respectively called refA, refB, refC and refD.

More information concerning the MAQC initiative can be found in the September 2006 special issue of *Nature Biotechnology*.

2 Loading the reference data

Once the library has been installed and loaded, the reference datasets can be loaded using the (data()) function as shown below.

```
> library("MAQCsubsetILM")
> data(refA)
> refA
Summary of data information:
Major Operation History:
            submitted
                                  finished
1 2008-02-29 12:24:41 2008-02-29 12:24:43
2 2008-02-29 12:24:43 2008-02-29 12:24:43
                                             command lumiVersion
1
     lumiR("ILM_1_A1.txt", parseColumnName = FALSE)
                                                           1.5.17
2 lumiQ(x.lumi = x.lumi, detectionTh = detectionTh)
                                                           1.5.17
. . .
             submitted
                                                                       command
                                   finished
72 2008-02-29 12:25:22 2008-02-29 12:25:22 combine(x = x.lumi, y = x.lumi.i)
73 2008-02-29 12:27:25 2008-02-29 12:27:25
                                                         Subsetting 3 samples.
   lumiVersion
72
        1.5.17
73
        1.5.17
Object Information:
LumiBatch (storageMode: lockedEnvironment)
assayData: 47293 features, 3 samples
  element names: beadNum, detection, exprs, se.exprs
```

```
protocolData: none
phenoData
  sampleNames: ILM_1_A5 ILM_2_A1 ILM_3_A2
  varLabels: sampleID site ref replicate
  varMetadata: labelDescription
featureData
  featureNames: GI_10047089-S GI_10047091-S ... trpF (47293 total)
  fvarLabels: TargetID
  fvarMetadata: labelDescription
experimentData: use 'experimentData(object)'
Annotation:
Control Data: N/A
QC information: Please run summary(x, 'QC') for details!
```