Package 'EpipwR'

January 2, 2025

Title Efficient Power Analysis for EWAS with Continuous or Binary

Description A quasi-simulation based approach to performing power analysis

Type Package

Version 1.0.0

Outcomes

```
for EWAS (Epigenome-wide association studies) with continuous or binary
      outcomes. 'EpipwR' relies on empirical EWAS datasets to determine power
      at specific sample sizes while keeping computational cost low. EpipwR can
      be run with a variety of standard statistical tests, controlling for either
      a false discovery rate or a family-wise type I error rate.
License Artistic-2.0
Encoding UTF-8
LazyData false
URL https://github.com/jbarth216/EpipwR
BugReports https://github.com/jbarth216/EpipwR
Imports EpipwR.data, ExperimentHub (>= 2.10.0), ggplot2
Depends R (>= 4.4.0)
Roxygen list(markdown = TRUE)
RoxygenNote 7.3.1
Suggests knitr, rmarkdown, testthat (>= 3.0.0), sessioninfo
VignetteBuilder knitr
biocViews Epigenetics, ExperimentalDesign
Config/testthat/edition 3
git_url https://git.bioconductor.org/packages/EpipwR
git_branch RELEASE_3_20
git_last_commit 9e2b25e
git_last_commit_date 2024-10-29
Repository Bioconductor 3.20
Date/Publication 2025-01-02
Author Jackson Barth [aut, cre] (<a href="https://orcid.org/0009-0009-6307-9928">https://orcid.org/0009-0009-6307-9928</a>),
      Austin Reynolds [aut],
      Mary Lauren Benton [ctb],
      Carissa Fong [ctb]
Maintainer Jackson Barth < Jackson_Barth@Baylor.edu>
```

2 EpipwR

Contents

	EpipwR_plot get_power_cc																		
	get_power_cont .	 			 •	•	•		•		 •			•	 •	•	•	•	5
Index																			7

Description

EpipwR: Efficient Power Analysis for EWAS with Continuous or Binary Outcomes

nary Outcomes

Details

A quasi-simulation based approach to performing power analysis for EWAS (Epigenome-wide association studies) with continuous or binary outcomes. 'EpipwR' relies on empirical EWAS datasets to determine power at specific sample sizes while keeping computational cost low. EpipwR can be run with a variety of standard statistical tests, controlling for either a false discovery rate or a family-wise type I error rate.

For a detailed overview, please use vignette("EpipwR").

Author(s)

Maintainer: Jackson Barth < Jackson_Barth@Baylor.edu > (ORCID)

Authors:

· Austin Reynolds

Other contributors:

- Mary Lauren Benton [contributor]
- Carissa Fong [contributor]

See Also

Useful links:

- https://github.com/jbarth216/EpipwR
- Report bugs at https://github.com/jbarth216/EpipwR

EpipwR_plot 3

EpipwR_plot

Plotting EpipwR Output

Description

Plots 95% error bars of EpipwR output. Different sample sizes are displayed along the x-axis with lines representing distinct correlation means. This functions should only be used with EpipwR output.

Usage

```
EpipwR_plot(df)
```

Arguments

df

A data frame containing output from the get_power_cont() or get_power_cc() function.

Value

A line plot with sample size on the x-axis and power on the y-axis. Multiple lines are plotted for each correlation/effect size mean specified.

Examples

```
out <- get_power_cont(100,10000,c(70,80,90,100),.05,c(.3,.4,.5)) 
 EpipwR_plot(out)
```

get_power_cc

Power Calculations for Case/Control EWAS

Description

Calculates power for EWAS with a binary outcome for multiple sample sizes and/or effect sizes. Data sets are only simulated for the non-null tests; p-values are generated directly for the null tests. Rather than specifying the number of data sets to calculate power, you specify the precision level (MOE) and a maximum number of data sets (Nmax). After 20 data sets, the function terminates when the desired precision level is reached or if the number of tested data sets reaches Nmax.

Usage

```
get_power_cc(
  dm,
  Total,
  n,
  fdr_fwer,
  delta_mu,
  delta_sd = 0,
  n1_prop = 0.5,
  Tissue = "Saliva",
```

4 get_power_cc

```
Nmax = 1000,
MOE = 0.03,
test = "pooled",
use_fdr = TRUE,
Suppress_updates = FALSE)
```

Arguments

Number of non-null tests.
The total number of tests (null and non-null).
Sample size(s) for which power is calculated (accepts a vector). This is total sample size.
Either the false discovery rate or the family-wise type I error rate, depending on use_fdr.
Average effect size (difference in mean methylation between groups) for non-null tests (accepts a vector). If multiple n and delta_mu are specified, power is calculated for all unique settings under the Cartesian product of these vectors.
Standard deviation of the effect size for non-null tests. If 0, all effect sizes are fixed at delta_mu.
Indicates the proportion of the total sample size (n) in group 1 (rounded to the nearest integer).
Tissue type of Empirical EWAS to be used for data generation (see details for valid options).
The maximum number of data sets used to calculate power.
The target margin of error of a 95% confidence interval for average power. This determines the stopping point of the algorithm.
The type of t-test to be used. "pooled" indicates a pooled variance t-test while "WS" indicates Welch's t-test.
If TRUE, uses fdr_fwer as the false discovery rate. If FALSE, uses the family-wise type I error rate.

Details

Suppress_updates

Valid options for the Tissue argument are "Saliva", "Lymphoma", "Placenta", "Liver", "Colon", "Blood adult", "Blood 5 year olds", "Blood newborns", "Cord-blood (whole blood)", "Cord-blood (PBMC)", "Adult (PBMC)", and "Sperm". All data sets are publicly available on the gene expression omnibus (see EpipwR.data package for more details) and were identified by Graw, et. al. (2019). Please note that, due to some extreme values in this data set, the Lymphoma option will occasionally throw a warning related to data generation. At this time, we recommend using one of the other tissue options.

If TRUE, blocks messages reporting the completion of each unique setting.

Unlike pwrEWAS (Graw et al., 2019), EpipwR enforces equality of precision (sum of the parameters) on the distributions of each group rather than equality of variance.

Value

A dataframe with rows equal to the number of n and rho_mu combinations

get_power_cont 5

References

Graw, S., Henn, R., Thompson, J. A., and Koestler, D. C. (2019). pwrEWAS: A user-friendly tool for comprehensive power estimation for epigenome wide association studies (EWAS). *BMC Bioinformatics*, 20(1):218.

Examples

```
# This examples calculates power for 100 non-null tests out of 10,000 total # with an FDR of 5%. Sample sizes of 40, 50, 60, 70 and fixed effect sizes # of 0.01, 0.02, 0.05 are used. For improved accuracy, MOE=.01 to ensure # that the 95% confidence interval for average power has a margin of error # no larger than .01 (unless the maximum of 1,000 data sets is exceeded). get_power_cc(100,10000,c(40,50,60,70),.05,c(.01,.02,.05), MOE=.01)
```

get_power_cont

Power Calculations for Continuous EWAS

Description

Calculates power for EWAS with a continuous outcome for multiple sample sizes and/or correlations. Data sets are only simulated for the non-null tests; p-values are generated directly for the null tests. Rather than specifying the number of data sets to calculate power, you specify the precision level (MOE) and a maximum number of data sets (Nmax). After 20 data sets, the function terminates when the desired precision level is reached or if the number of tested data sets reaches Nmax.

Usage

```
get_power_cont(
   dm,
   Total,
   n,
   fdr_fwer,
   rho_mu,
   rho_sd = 0,
   Tissue = "Saliva",
   Nmax = 1000,
   MOE = 0.03,
   test = "pearson",
   use_fdr = TRUE,
   Suppress_updates = FALSE
)
```

Arguments

dm Number of non-null tests.

Total The total number of tests (null and non-null).

n Sample size(s) for which power is calculated (accepts a vector).

fdr_fwer Either the false discovery rate or the family-wise type I error rate, depending on

 $use_fdr.$

6 get_power_cont

rho_mu	Mean correlation(s) of methylation and phenotype for non-null tests (accepts a vector). If multiple n and rho_mu are specified, power is calculated for all unique settings under the Cartesian product of these vectors.
rho_sd	Standard deviation of methylation and phenotype for non-null tests. If 0 all correlations are fixed at rho_mu.
Tissue	Tissue type of Empirical EWAS to be used for data generation (see details for valid options).
Nmax	The maximum number of datasets used to calculate power
MOE	The target margin of error of a 95% confidence interval for average power. This determines the stopping point of the algorithm.
test	The type of statistical test to be used. "pearson", "kendall" and "spearman" are valid options.
use_fdr	If TRUE, uses fdr_fwer as the false discovery rate. If FALSE, uses the family-wise type I error rate.

Suppress_updates

If TRUE, blocks messages reporting the completion of each unique setting.

Details

Valid options for the Tissue argument are "Saliva", "Lymphoma", "Placenta", "Liver", "Colon", "Blood adult", "Blood 5 year olds", "Blood newborns", "Cord-blood (whole blood)", "Cord-blood (PBMC)", "Adult (PBMC)", and "Sperm". All data sets are publicly available on the gene expression omnibus (see EpipwR.data package for more details) and were identified by Graw, et. al. (2019). Please note that, due to some extreme values in this data, the Lymphoma option will occasionally throw a warning related to data generation. At this time, we recommend using one of the other tissue options.

Although this function only covers 3 types of statistical tests, its worth noting that tests run using software packages such as limma will yield the same results as a pearson correlation test in the absence of covariates or any dependence across CpG sites (as is the assumption here). Any users wanting to mimic an analysis done in limma should use test="pearson".

Value

A data frame with rows equal to the number of n and rho_mu combinations

References

Graw, S., Henn, R., Thompson, J. A., and Koestler, D. C. (2019). pwrEWAS: #' A user-friendly tool for comprehensive power estimation for epigenome wide association studies (EWAS). *BMC Bioinformatics*, 20(1):218.

Examples

```
# This examples calculates power for 100 non-null tests out of 10,000 total # with an FDR of 5%. # Sample sizes of 70,80,90,100 and fixed correlations at 0.3,0.4,0.5 are # used. get_power_cont(100,10000,c(70,80,90,100),.05,c(.3,.4,.5))
```

Index