Package 'flowType'

October 12, 2016

Type Package

| Title Phenotyping Flow Cytometry Assays |
|---|
| Version 2.10.0 |
| Date 2014-11-26 |
| Author Nima Aghaeepour, Kieran O'Neill, Adrin Jalali |
| Maintainer Nima Aghaeepour <naghaeep@gmail.com></naghaeep@gmail.com> |
| Description Phenotyping Flow Cytometry Assays using multidimentional expansion of single dimentional partitions. |
| Imports Biobase, graphics, grDevices, methods, flowCore, flowMeans, sfsmisc, rrcov, flowClust, flowMerge, stats |
| Depends R (>= 2.10), Rcpp (>= 0.10.4), BH (>= 1.51.0-3) |
| LinkingTo Rcpp, BH |
| Suggests xtable |
| biocViews FlowCytometry |
| License Artistic-2.0 |
| LazyLoad yes |
| NeedsCompilation yes |
| R topics documented: |
| flowType-package |
| calcMemUse 3 |
| calcNumPops |
| decodePhenotype |
| DLBCLExample |
| encodePhenotype 6 flowType 7 |
| getLabels |
| HIVData |
| HIVMetaData |
| Phenotypes-class |
| plot |
| summary-methods |

2 flowType-package

Index 15

flowType-package flowType: Phenotyping Flow Cytometry Assays

Description

flowType uses a simple threshold, Kmeans, flowMeans, or flowClust to partition every channel to a positive and a negative cell population. These partitions are then combined to generate a set of multi-dimensional phenotypes.

Details

Package: flowType
Type: Package
Version: 0.0.1
Date: 2011-04-27
License: Artistic-2.0
LazyLoad: yes
Depends: methods

For a given FCS file, the flowType function extracts a the phenotypes and reports their cell frequencies (number of cells) and mean fluorescence intensity (MFI)s.

Author(s)

Nima Aghaeepour, Kieran O'Neill, Adrin Jalali

References

Please cite the following for the current version of flowType:

O'Neill K, Jalali A, Aghaeepour N, Hoos H, Brinkman RR. Enhanced flowType/RchyOptimyx: a BioConductor pipeline for discovery in high-dimensional cytometry data. Bioinformatics. 2014 May 1;30(9):1329-30. doi: 10.1093/bioinformatics/btt770

The original paper and description can be found at:

Nima Aghaeepour, Pratip K. Chattopadhyay, Anuradha Ganesan, Kieran O'Neill, Habil Zare, Adrin Jalali, Holger H. Hoos, Mario Roederer, and Ryan R. Brinkman. Early Immunologic Correlates of HIV Protection can be Identified from Computational Analysis of Complex Multivariate T-cell Flow Cytometry Assays. Bioinformatics, 2011.

```
#Load the library
library(flowType)
data(DLBCLExample)
MarkerNames <- c('Time', 'FSC-A','FSC-H','SSC-A','IgG','CD38','CD19','CD3','CD27','CD20', 'NA', 'NA')</pre>
```

calcMemUse 3

```
#These markers will be analyzed
PropMarkers <- 3:5
MFIMarkers <- PropMarkers
MarkerNames <- c('FS', 'SS', 'CD3', 'CD5', 'CD19')</pre>
#Run flowType
Res <- flowType(DLBCLExample, PropMarkers, MFIMarkers, 'kmeans', MarkerNames);</pre>
MFIs=Res@MFIs;
Proportions=Res@CellFreqs;
Proportions <- Proportions / max(Proportions)</pre>
names(Proportions) <- unlist(lapply(Res@PhenoCodes,</pre>
                       function(x){return(decodePhenotype(
                       x,Res@MarkerNames[PropMarkers],
                       Res@PartitionsPerMarker))}))
#Select the 30 largest phenotypes
index=order(Proportions, decreasing=TRUE)[1:30]
bp=barplot(Proportions[index], axes=FALSE, names.arg=FALSE)
text(bp+0.2, par("usr")[3]+0.02, srt = 90, adj = 0, labels = names(Proportions[index]), xpd = TRUE, cex=0.8)
axis(2);
axis(1, at=bp, labels=FALSE);
title(xlab='Phenotype Names', ylab='Cell Proportion')
#These phenotype can be analyzed using a predictive model (e.g., classification or regression)
```

calcMemUse

Function calcMemUse in Package flowType

Description

Estimates the memory usage in bytes for running flowType with a given set of parameters.

Usage

calcMemUse(NumPops, NumPropMarkers, NumMFIMarkers, NumCells, MaxMarkersPerPop, PartitionsPerChannel

Arguments

NumPops Number of cell types which will be returned. Can be computed using calcNumPops

NumPropMarkers Numer of markers to use for combinatorial gating

NumMFIMarkers Number of markers to determine the MFIs of for every cell type

NumCells Number of cells in the flowFrame passed to flowType

MaxMarkersPerPop

Maximum number of markers to use at once in combinatorial gating (ie all cell

types over 1:MaxMarkersPerPop will be counted)

PartitionsPerChannel

Number of partitions per channel.

4 calcNumPops

Details

If you use different numbers of partitions for different channels, try providing the highest number as PartitionsPerChannel, and expect an over-estimate.

Value

Estimated memory use in bytes.

Author(s)

Kieran O'Neill

See Also

calcNumPops, flowType

calcNumPops

Function calcNumPops in Package flowType

Description

Compute the number of populations that will be produced by running flowType with a given set of parameters. This is especially useful for estimating memory but is also used internally to determine the size of return objects to pass down to C++. This may also be useful for determining the cutoff for number of markers to use to make phenotypes (in terms of statistical power for later testing).

Usage

```
calcNumPops(PartitionsPerMarker, MaxMarkersPerPop)
```

Arguments

PartitionsPerMarker

Integer vector specifying the number of partitions for each marker, in order.

MaxMarkersPerPop

Integer speciying the threshold chosen

Value

Integer specifying the number of populations the given paramters would produce.

Author(s)

Kieran O'Neill

See Also

calcMemUse, flowType

decodePhenotype 5

Examples

```
calcNumPops(c(2,2,3,2,2,4), 5)
```

decodePhenotype

Method decodePhenotype *in Package* **flowType**

Description

Method to decode phenotypes back to a human-readable string.

Details

FlowType's encoding is as follows:

0 – marker not considered in phenotype (don't care about its value) 1 – marker is negative (e.g. CD4-) 2 – marker is positive (e.g. CD4+) 3 – marker is positive, but brighter than 2 (CD4++) 4 – marker is even brighter (CD4+++) etc

Note that this encoding system does not allow for "dim" markers – dim positivity is denoted by the first level of positivity.

Also note that the encoding is performed from the dimmest to the brightest partition, but the location of thresholds will dictate the interpretation of the code. (e.g. if you only set one threshold, but you place it between the positive and the bright population, then both positive and negative events will be considered negative.)

Methods

```
signature(pheno.code = "character", marker.names = "character", partitions.per.marker = "numeric")
In flowType, phenotypes themselves are represented by codes (e.g. 012) rather than full
strings (CD4+CD8-), in order to save memory when a very large number of phenotypes are
considered. decodePhenotype serves to translate the codes back to a human-readable string.
```

Author(s)

Kieran O'Neill

See Also

```
encodePhenotype, flowType
```

```
decodePhenotype('1034',c('CD34','CD3','CD45','CD19'), 4)
decodePhenotype('20013',c('CD34','CD3','CD45','CD19', 'CD20'), 4)
```

6 encodePhenotype

DLBCLExample

DLBCLExample

Description

A flow cytometry sample from a patient with DLBC lymphoma. The full dataset is available through the FlowCAP project (http://flowcap.flowsite.org).

Usage

```
data(DLBCLExample)
```

Format

A flowFrame describing expression values of 3 markers and 3796 cells. Each column represents a marker and each row represents a cell.

Examples

data(DLBCLExample)

encodePhenotype

Function encodePhenotype in Package flowType

Description

In **flowType**, phenotypes themselves are represented by codes (e.g. 012) rather than full strings (CD4+CD8-), in order to save memory when a very large number of phenotypes are considered. encodePhenotype serves to translate a human-readable string down to flowType's internal coded representation.

Usage

```
encodePhenotype(pheno.string, marker.names)
```

Arguments

pheno.string character vector containing containing the phenotype string to be encoded marker.names vector of character vectors each specifying the name of a channel, in order

flowType 7

Details

FlowType's encoding is as follows:

0 – marker not considered in phenotype (don't care about its value) 1 – marker is negative (e.g. CD4-) 2 – marker is positive (e.g. CD4+) 3 – marker is positive, but brighter than 2 (CD4++) 4 – marker is even brighter (CD4+++) etc

Note that this encoding system does not allow for "dim" markers – dim positivity is denoted by the first level of positivity.

Also note that the encoding is performed from the dimmest to the brightest partition, but the location of thresholds will dictate the interpretation of the code. (e.g. if you only set one threshold, but you place it between the positive and the bright population, then both positive and negative events will be considered negative.)

Value

Character vector containing the encoded phenotype, with one character per channel.

Author(s)

Kieran O'Neill

See Also

decodePhenotype, flowType

Examples

```
encodePhenotype('CD34++CD3-CD45+++',c('CD34','CD3','CD45','CD19'))
```

flowType

flowType: Phenotyping Flow Cytometry Assays

Description

flowType uses a simple threshold, Kmeans, flowMeans or flowClust to partition every channel to a positive and a negative cell population. These partitions are then combined to generate a set of multi-dimensional phenotypes.

Usage

flowType(Frame, PropMarkers=NULL, MFIMarkers=NULL, Methods='kmeans', MarkerNames=NULL, MaxMarkersPer #If upgrading from flowType 1.x to 2.x, please check documentation as some arguments have changed slig

8 flowType

Arguments

Frame A flowFrame (after transformation) that is going to be phenotyped.

PropMarkers A vector of the indexes or names of the markers to partition to specify pheno-

types. If NULL, all markers in the frame will be used.

MFIMarkers A vector of the indexes or names of the markers for which MFIs must be mea-

sured. If NULL, no markers will be used.

Methods A single string specifying the method to use to determine thresholds for par-

titioning of markers. Values can be "kmeans", "flowMeans", "flowClust", or "Thresholds". If "Thresholds" is specified, user-specified thresholds must be

provided via the Thresholds parameter.

MarkerNames A vector of names for the channels. If NULL, the names in Frame will be used.

MaxMarkersPerPop

An integer specifying the maximum number of markers to use to define populations (how "deep" to phenotype). This should be less than or equal to PropMarkers.

If NULL, will default to the length of PropMarkers.

PartitionsPerMarker

An integer or vector of integers specifing the number of partitions per marker. If a single integer, this number will be used for all markers. If a vector, the

numbers will be matched with PropMarkers in order.

Thresholds A list of vectors specifying per-channel thresholds. Each list item corresponds

to one marker, and contains the threshold(s) for that marker. If only one vector is provided in the list, then those thresholds will be used for all markers. Otherwise, the list must be of the same length as PropMarkers. Note: if Methods == 'thresholds', then Thresholds must be specified. If not, it is

ignored.

MemLimit Memory limit in GB. flowType will do a sanity check before executing, and

if the total size of counts plus MFI values for all populations would exceed

MemLimit, will not run.

verbose Boolean variable. If TRUE, information about different processing tasks will be

printed into the standard output.

Value

CellFreqs: Object of class "numeric" containing the cell frequencies measured for each

phenotype. Phenotype names are assigned as labels.

MFIs: Object of class "matrix" containing the measured MFIs for each phenotype.

Phenotype names are assigned as column labels and marker names as row labels.

PropMarkers A vector of the indexes or names of the markers for which cell proportions must

be measured.

MFIMarkers A vector of the indexes or names of the markers for which MFIs must be mea-

sured.

MarkerNames A vector of names for the channels. If NULL, the names provided in Frame will

be used.

Partitions A matrix where each column shows the partitioning of the respective channel. 1

and 2 correspond to negative and positive, respectively.

flowType 9

PhenoCodes

A vector of strings of length N (the number of markers) for each phenotype measured. For every phenotype, the character corresponding to a given marker can be 0, 1, 2, etc for neutral, negative, positive, bright, etc. See the provided vigenette for more details and examples.

Author(s)

Nima Aghaeepour, Kieran O'Neill

References

Please cite the following for the current version of flowType:

O'Neill K, Jalali A, Aghaeepour N, Hoos H, Brinkman RR. Enhanced flowType/RchyOptimyx: a BioConductor pipeline for discovery in high-dimensional cytometry data. Bioinformatics. 2014 May 1;30(9):1329-30. doi: 10.1093/bioinformatics/btt770

The original paper and description can be found at:

Nima Aghaeepour, Pratip K. Chattopadhyay, Anuradha Ganesan, Kieran O'Neill, Habil Zare, Adrin Jalali, Holger H. Hoos, Mario Roederer, and Ryan R. Brinkman. Early Immunologic Correlates of HIV Protection can be Identified from Computational Analysis of Complex Multivariate T-cell Flow Cytometry Assays. Bioinformatics, 2011.

```
#Load the library
library(flowType)
data(DLBCLExample)
#These markers will be analyzed
PropMarkers <- 3:5
MFIMarkers <- PropMarkers
MarkerNames <- c('FS', 'SS', 'CD3', 'CD5', 'CD19')</pre>
#Run flowType
Res <- flowType(DLBCLExample, PropMarkers, MFIMarkers, 'kmeans', MarkerNames);</pre>
MFIs=Res@MFIs;
Proportions=Res@CellFreqs;
Proportions <- Proportions / max(Proportions)</pre>
names(Proportions) <- unlist(lapply(Res@PhenoCodes,</pre>
                       function(x){return(decodePhenotype(
                       x, Res@MarkerNames[PropMarkers],
                       Res@PartitionsPerMarker))}))
#Select the 30 largest phenotypes
index=order(Proportions, decreasing=TRUE)[1:30]
bp=barplot(Proportions[index], axes=FALSE, names.arg=FALSE)
text(bp+0.2, par("usr")[3]+0.02, srt = 90, adj = 0, labels = names(Proportions[index]), xpd = TRUE, cex=0.8)
axis(2);
axis(1, at=bp, labels=FALSE);
title(xlab='Phenotype Names', ylab='Cell Proportion')
```

10 getLabels

#These phenotype can be analyzed using a predictive model (e.g., classification or regression)

getLabels

getLabels: Returns the labels of the cells in a given phenotype.

Description

Returns the labels of the cells in a given phenotype in a Phenotypes object.

Usage

```
getLabels(Phenotypes, PhenotypeNumber)
```

Arguments

Phenotypes

An object of class Phenotypes as produced by the flowType function.

PhenotypeNumber

A numeric or character value representing the phenotypes number of name, re-

spectively.

Value

Membership Labels:

A vector of length of the number of events. 1 and 2 represent the cells that are not-included and included in the phenotype respectively.

Author(s)

Nima Aghaeepour

References

Nima Aghaeepour, Pratip K. Chattopadhyay, Anuradha Ganesan, Kieran O'Neill, Habil Zare, Adrin Jalali, Holger H. Hoos, Mario Roederer, and Ryan R. Brinkman. Early Immunologic Correlates of HIV Protection can be Identified from Computational Analysis of Complex Multivariate T-cell Flow Cytometry Assays. submitted to Bioinformatics, 2011.

Examples

#See the vigentte

HIVData 11

HIVData HIVData

Description

A flow cytometry dataset from a HIV+ patients PBMC by the Scott lab of the Simon Fraser University and the Spina lab of the University of California San Diego.

Usage

data(HIVData)

Format

A flowSet describing expression values of 11 markers and 500 cells (sampled uniformly) for 19 HIV+ and 12 normal subjects.

Examples

data(HIVData)

HIVMetaData

HIVMetaData

Description

The meta-data of a flow cytometry dataset from a HIV+ patients PBMC by the Scott lab of the Simon Fraser University and the Spina lab of the University of California San Diego.

Usage

data(HIVMetaData)

Format

A matrix describing the FCS filename, patient label (HIV+ or normal) and tube number of every assay.

Examples

data(HIVMetaData)

Phenotypes-class

Phenotypes-class

Class "Phenotypes"

Description

The return data from running flowType, containing counts

Objects from the Class

Objects can be created by calls of the form new("Phenotypes", ...).

Slots

CellFreqs: Numeric vector containing counts of the number of cells belonging to each cell type.

MFIs: Matrix of MFIs, with rows for cell types and columns for markers.

PhenoCodes: Vector of character strings representing the codes of each cell type (phenotype).

PropMarkers: Numeric vector specifying which markers were used for combinatorial gating.

MFIMarkers: Numeric vector specifying for which markers MFIs were computed for each cell type.

MarkerNames: A character vector of the names of all markers in the flowFrame given

Partitions: The first level partitions that each cell in the flowFrame belong to in each channel.

MaxPopSize: MaxMarkersPerPop

PartitionsPerMarker: Vector of number of partitions used for each marker

Thresholds: A list of vectors with the calculated thresholds for each marker (if a clustering algorithm was used) or the thresholds provided by the user.

See Also

flowType

```
showClass("Phenotypes")
```

plot 13

plot

Methods for Function plot

Description

Methods for function plot

Usage

```
## S4 method for signature 'Phenotypes,flowFrame'
plot(x, y, ...)
## S4 method for signature 'Phenotypes,numeric'
plot(x, y, Frame,...)
## S4 method for signature 'Phenotypes,character'
plot(x, y, Frame,...)
```

Arguments

| X | An object of class Phenotypes as generated by the flowType package. |
|-------|---|
| у | A flowFrame or a numeric/character value representing the phenotype that needs to be plotted depending on the signature of the function |
| Frame | A flowFrame (might be optional depending on the signature of the function |
| | Extra parameters that will be passed to the generic plot function |

Author(s)

Nima Aghaeepour <<naghaeep@gmail.com>>

See Also

flowType

```
#See the vigentte
```

summary-methods

 $\verb"summary-methods"$

~~ Methods for Function summary ~~

Description

 $\sim\sim$ Methods for function summary $\sim\sim$

Methods

signature(object = "Phenotypes") Prints basic characteristics of a Phenotypes object.

See Also

flowType

Examples

#See the vigentte

Index

| *Topic FlowCytData | decodePhenotype-methods |
|---|--|
| flowType, 7 | (decodePhenotype), 5 |
| flowType-package, 2 | DLBCLExample, 6 |
| getLabels, 10 | |
| *Topic HIV | encodePhenotype, 5, 6 |
| flowType, 7 | encodePhenotype,character,character-method |
| flowType-package, 2 | (encodePhenotype), 6 |
| getLabels, 10 | encodePhenotype-methods |
| *Topic classes | (encodePhenotype), 6 |
| Phenotypes-class, 12 | |
| *Topic classification | flowType, 4, 5, 7, 7, 12–14 |
| flowType, 7 | flowType-package, 2 |
| flowType-package, 2 | gott shallo 10 |
| getLabels, 10 | getLabels, 10 |
| *Topic clustering | HIVData, 11 |
| flowType, 7 | HIVMetaData, 11 |
| flowType-package, 2 | Titric cabaca, 11 |
| getLabels, 10 | Phenotypes (Phenotypes-class), 12 |
| *Topic datasets | Phenotypes-class, 12 |
| DLBCLExample, 6 | plot, 13 |
| HIVData, 11 | plot, Phenotypes, character (plot), 13 |
| HIVMetaData, 11 | plot,Phenotypes,character-method |
| *Topic print | (plot), 13 |
| plot, 13 | plot, Phenotypes, flowFrame (plot), 13 |
| summary-methods, 14 | plot, Phenotypes, flowFrame-method |
| *Topic utilities | (plot), 13 |
| calcMemUse, 3 | plot, Phenotypes, numeric (plot), 13 |
| calcNumPops, 4 | <pre>plot,Phenotypes,numeric-method(plot),</pre> |
| decodePhenotype, 5 | 13 |
| encodePhenotype, 6 | |
| | summary,Phenotypes-method |
| calcMemUse, 3, 4 | (summary-methods), 14 |
| calcNumPops, 3, 4, 4 | summary-methods, 14 |
| decodePhenotype, 5, 7 | |
| decodePhenotype,character, | |
| character, numeric-method | |
| (decodePhenotype), 5 | |
| decodePhenotype,character,character,num | eric-method |
| (decodePhenotype), 5 | |