Package 'flowQ'

October 17, 2017

Type Package

Title Quality control for flow cytometry
Version 1.36.0
Author R. Gentleman, F. Hahne, J. Kettman, N. Le Meur, N. Gopalakrishnan
Maintainer Mike Jiang <wjiang2@fhcrc.org></wjiang2@fhcrc.org>
Description Provides quality control and quality assessment tools for flow cytometry data.
Depends R (>= 2.10.0), methods, BiocGenerics, outliers, lattice, flowViz, myoutlier, bioDist, parody, RColorBrewer, latticeExtra
Imports methods, BiocGenerics, geneplotter, flowCore, flowViz, IRanges
Suggests flowStats
SystemRequirements ImageMagick
License Artistic-2.0
biocViews Infrastructure, FlowCytometry, CellBasedAssays
NeedsCompilation no
R topics documented:
flowQ-package
aggregatorList-class
binaryAggregator-class
discreteAggregator-class
evaluateProcess
factorAggregator-class

 locateDuplicatedParameters
 8

 normailzeSets
 8

 numericAggregator-class
 9

 outlier-class
 10

 outlierResult-class
 11

 outliers-methods
 11

 preProcessFlowList
 12

 qaAggregator-class
 12

 qaGraph-class
 13

 qaGraphList-class
 14

 qaProcess-class
 15

 qaProcess.2DStatsPlot
 16

5 6 2 flowQ-package

flow	Quality control for flow cytometry																			
Index																				38
	writeQATextReport			• •				•			•			•	 •	 •	•	•	•	. 30
	writeQAReport																			
	validProcess																			
	txtFormatQAObject																			
	stringAggregator-cl																			
	rangeAggregator-cla	ass																		. 31
	qData																			
	qaReport																			
	qaProcessSummary																			
	qaProcessFrame-cla																			
	qaProcess.timeline																			
	qaProcess.timeflow																			
	qaProcess.marginev																			
	qaProcess.KLDistPl																			
	qaProcess.DensityP qaProcess.ECDFPlo																			
	qaProcess.cellnumb																			
	qaProcess.Boundary																			

Description

Functions and methods for quality control and QA of flow cytometry data. This package heavily depends on the flowCore package.

Details

Package: flowQ Type: Package Version: 0.2.1

Date: 2006-11-16 License: Artistic

Quality control is an important aspect when dealing with large amounts of complex high-throughput data. This package comprises functionality for efficient QA of flow cytometry data.

Author(s)

Maintainer: Florian Hahne <f.hahne@dkfz.de> Authors: R. Gentleman, F. Hahne, J. Kettman, N. Le Meur, M. Tang

References

references go here

aggregatorList-class 3

See Also

flowCore

Examples

```
## examples go here
```

```
aggregatorList-class Class "aggregatorList"
```

Description

A list of qaAggregators

Details

This class directly extends class "list" and is intended to exclusively hold objects of class qaAggregator, where each list item represents the outcome of a QA subprocess for a single flowFrame. It mainly exists to allow for method dispatch and should never be populated manualy; instead, use the constructor qaAggregatorList which checks for valid objects.

Objects from the Class

Objects should be created using the constructor:

```
aggregatorList(...), where ... are objects inheriting form class qaAggregator or a list of such objects.
```

Slots

```
.Data: Object of class "list", the list data
```

Extends

```
Class "list", from data part. Class "vector", by class "list", distance 2.
```

Methods

```
initialize signature(.Object = "aggregatorList"): constructor
show signature(object = "aggregatorList"): print object details
```

Author(s)

Florian Hahne

See Also

```
qaGraph, writeQAReport, qaProcess, qaAggregator
```

```
showClass("aggregatorList")
```

binaryAggregator-class

Class "binaryAggregator"

Description

Abstraction of a binary type of aggregator with possible states "passed" and "not passed"

Objects from the Class

Objects can be created by calls of the form new("binaryAggregator", ...), or using the constructor binaryAggregator(passed), where passed is a logical scalar.

Slots

passed: Object of class "logical" indicating whether the process has passed the QA requirements

Extends

Class qaAggregator, directly.

Methods

```
show signature(object = "binaryAggregator"): print object details
writeLines signature(text = "binaryAggregator", con = "file", sep = "missing", useBytes="missing"
write to HTML file connection
```

Author(s)

Florian Hahne

See Also

```
qaProcess.marginevents, qaReport, qaProcess, qaProcess.timeline, discreteAggregator, factorAggregator, numericAggregator, stringAggregator, rangeAggregator
```

```
showClass("binaryAggregator")
```

discreteAggregator-class

5

```
discreteAggregator-class
```

Class "discreteAggregator"

Description

Abstraction of a discrete type of aggregator with possible states "passed", "not passed" and "warning"

Objects from the Class

Objects can be created by calls of the form new("discreteAggregator", ...) or using the constructor discreteAggregator(x), where x is an integer value in [0,1,2] or a factor with levels 0, 1 and 2.

Slots

```
x: Object of class "factor" One in 0 (not passes), 1 (passed) or 2 (warning)
passed: Object of class "logical" indicating whether the process has passed the QA requirements
(not eveluated)
```

Extends

Class qaAggregator, directly.

Methods

Author(s)

Florian Hahne

See Also

```
qaProcess.marginevents, qaReport, qaProcess, qaProcess.timeline, binaryAggregator, factorAggregator, numericAggregator, stringAggregator, rangeAggregator
```

```
showClass("discreteAggregator")
```

6 evaluateProcess

evaluateProcess

Evaluate QA processes

Description

Re-evaluate an object of class qaProcess, e.g. for the case that a threshold value has changed.

Usage

```
evaluateProcess(process, thresh, ...)
```

Arguments

process An object of class qaProcess.

thresh The new treshold on which the process is to be evaluated.

... Further arguments that are passed on to the individual functions for each QA

process type.

Details

It is sometimes useful to update the state of aggregators in a qaProcess, for instances after changing the threshold value, without having to recompute all images, which can be very time consuming.

Value

An updated object of class qaProcess

Note

This function needs to be extended for new types of qaProcess.

Author(s)

Florian Hahne

See Also

```
qaProcess, writeQAReport
```

factorAggregator-class 7

```
factorAggregator-class
```

Class "factorAggregator"

Description

Abstraction of a factor type of aggregator with possible states coded by the factor levels

Objects from the Class

Objects can be created by calls of the form new("factorAggregator", ...) or using the constructor factorAggregator(x, passed), where x is a factor, or an object which can be coerced to a factor, and passed is a logical scalar.

Slots

```
x: Object of class "factor" coding the outcome state
passed: Object of class "logical" indicating whether the process has passed the QA requirements
```

Extends

```
Class "qaAggregator", directly.
```

Methods

Author(s)

Florian Hahne

See Also

```
qaProcess.marginevents, qaReport, qaProcess, qaProcess.timeline, discreteAggregator, binaryAggregator, numericAggregator, stringAggregator, rangeAggregator,
```

```
showClass("factorAggregator")
```

8 normailzeSets

locateDuplicatedParameters

Locate Duplicated Parameters

Description

Identifies dyes that are repeated arcross multiple flowSets using information provided in the description field of each flowFrame.

Usage

locateDuplicatedParameters(flowList)

Arguments

flowList A list of flowSet

Details

This function is used internally by the qaProcess functions to identify dyes that are repeated across multiple aliquots. The function takes in a list of flowSet and identifies stains that are repeated across multiple flowSets (aliquots). These parameters are expected to have the same distribution across aliquots and can be used for Quality control. The method looks for the information provided in the description field of each flowFrame.

Value

An object of class qaProcess.

Author(s)

Nishant Gopalakrishnan

See Also

normalizeSets

normailzeSets

Locate Duplicated Parameters

Description

Aligns up peaks for parameters across multiple flowSets using the warpSet method from the flowStats package.

Usage

normalizeSets(flowList,dupes,peaks=NULL)

Arguments

flowList A list of flowSet

dupes Stains that are to be aligned across multiple flowSet

peaks Number of peaks expected in the distribution. If NULL is passed as input, The

number of peaks to align are estimated from the distribution

Details

This function can be used to normalize dye information across multiple flowSets so that the peaks in the distributions align. The method makes use of information provided in the description field of each flowFrame to align up floursecense nformation from multiple flowSet.

Value

A list of objects of class flowSet.

Author(s)

Nishant Gopalakrishnan

See Also

locateDuplicatedParameters warpSet

numericAggregator-class

Class "numericAggregator"

Description

Abstraction of a numeric type of aggregator for which possible states are coded by a numeric value

Objects from the Class

Objects can be created by calls of the form new("numericAggregator", ...) or using the constructor numericAggregator(x, passed), where x is a numeric scalar, and passed is a logical scalar.

Slots

x: Object of class "numeric" coding the outcome state.

passed: Object of class "logical" indicating whether the process has passed the QA requirements.

Extends

Class qaAggregator, directly.

10 outlier-class

Methods

Author(s)

Florian Hahne

See Also

qa Process. margine vents, qa Report, qa Process, qa Process. timeline, discrete Aggregator, factor Aggregator, binary Aggregator, string Aggregator, range Aggregator

Examples

```
showClass("numericAggregator")
```

outlier-class

Virtual Class "outlier"

Description

A class to represent outlier tests

Objects from the Class

A virtual Class: No objects may be created from it.

Slots

```
test: Object of class "character" ~~ parameters: Object of class "ANY" ~~
```

Methods

No methods defined with class "outlier" in the signature.

References

~put references to the literature/web site here ~

```
## showClass("outlier")
```

outlierResult-class 11

Description

A class to hold the results of an outlier test

Objects from the Class

```
Objects can be created by calls of the form new("outlierResult", ...). ~~ describe objects here ~~
```

Slots

```
frameId: Object of class "character" ~~
filterDetails: Object of class "list" ~~
test: Object of class "character" ~~
parameters: Object of class "ANY" ~~
```

Extends

```
Class "outlier", directly.
```

Methods

No methods defined with class "outlierResult" in the signature.

References

~put references to the literature/web site here ~

Examples

```
#showClass("outlierResult")
```

```
outliers-methods
```

~~ Methods for Function outliers in Package '.GlobalEnv' ~~

Description

```
{\sim}{\sim} Methods for function outliers in Package '.GlobalEnv' {\sim}{\sim}
```

Methods

```
x = "flowSet" Desribe this method here
```

12 qaAggregator-class

preProcessFlowList Preprocess a list of flowSets to add empty flowFrames for patients in which certain aliquots are missing

Description

The preProcessFlowList function takes in a list of flowSets as inputs and adds empty flowFrames for patients in which data from certain aliquots are missing. The function is to be used to modify datasets to conform to a format that can be utilized by the multiple sample qaProcess functions.

Usage

preProcessFlowList(flowList)

Arguments

flowList

A list of flowSet in which flowFrames for certain patients are missing

Details

This function can be used to modify a list of flowSets to fill in empty flowFrames for patients for which samples are missing.

Value

A list of objects of class flowSet.

Author(s)

Nishant Gopalakrishnan

See Also

locateDuplicatedParameters

qaAggregator-class

Abstraction of the possible outcomes of a QA process

Description

Virtual parent class for different types of QA aggregators

Details

In the context of this package, qaAggregators are objects that hold the outcome of a QA process. Each subclass implements its own writeLine method, which creates the appropriate HTML code for a graphical representation of the object.

Objects from the Class

A virtual Class: No objects may be created from it.

qaGraph-class 13

Slots

passed: Object of class "logical" indicating whether the process has passed the QA requirements

Methods

No methods defined with class "qaAggregator" in the signature.

Author(s)

Florian Hahne

See Also

qaProcess.marginevents, qaReport, qaProcess, qaProcess.timeline, binaryAggregator, discreteAggregator, factorAggregator, numericAggregator, stringAggregator, rangeAggregator

Examples

```
showClass("qaAggregator")
```

qaGraph-class

Class "qaGraph"

Description

Abstraction of the graphical output created in the cause of a QA process

Objects from the Class

Objects should be created using the constructor:

qaGraph(fileName, imageDir, width, empty=FALSE, pdf=TRUE), where fileName is a path to an image file, imageDir is the destination path for the images and the optional argument width is the final width to which the bitmap images are converted. The pdf argument controls wether vectorized versions of the image are to be produced. For the special case of an empty object without any images one can use option empty=TRUE, in which case the constructor ignores all other arguments.

During object instantiation the image file will be converted, resized and copied if necessary.

Slots

fileNames: Object of class "character" The paths to the image files, both the vectorized and unverctorized versions

dimensions: Object of class "matrix" The dimensions of the image files, both for the vectorized and unvectorized version

types: Object of class "character" The file extensions for both versions

id: Object of class "character" A unique identifier for the images

Methods

```
initialize signature(.Object = "qaGraph"): constructor
names signature(x = "qaGraph"): returns the file name of the bitmap version of the image
show signature(object = "qaGraph"): print object details
```

14 qaGraphList-class

Author(s)

Florian Hahne

See Also

```
qaGraphList, writeQAReport, qaProcess
```

Examples

```
showClass("qaGraph")
```

qaGraphList-class

Class "qaGraphList"

Description

A list of qaGraph objects

Details

This class directly extends class "list" and is intended to exclusively hold objects of class qaGraph, where each list item represents the grapical output of a QA subprocess for a single flowFrame. It mainly exists to allow for method dispatch and should never be populated manualy; instead, use the constructor qaGraphList which makes sure, that all image files are converted into the appropriate types and sizes and copied to the expected file location.

Objects from the Class

Objects should be created using the constructor:

qaGraphList(imageFiles, imageDir, width, pdf=TRUE), where imageFiles are paths to image files, imageDir is the destination path for the images and width is the final width to which the bitmap images are converted. pdf can be used to toggle the generation of vectorized versions of the images.

Slots

```
.Data: Object of class "list", the list data
```

Extends

```
Class "list", from data part. Class "vector", by class "list", distance 2.
```

Methods

```
initialize signature(.Object = "qaGraphList"): constructor
show signature(object = "qaGraphList"): print object details
```

Author(s)

Florian Hahne

qaProcess-class 15

See Also

```
qaGraph, writeQAReport, qaProcess
```

Examples

```
showClass("qaGraphList")
```

qaProcess-class

Abstraction of the results of a QA process

Description

QA processes create graphical output which can be bundled in a single HTML document. This class stores all information that is needed by writeQAReport to produce such HTML reports.

Objects from the Class

Objects should be created using the constructor functions. See qaProcess.timeline and qaProcess.marginevents for details. When writing new QA process functions, the constructors qaProcessFrame and qaProcess should be used. The latter expects the mandatory arguments id, type and frameProcesses and also accepts the optional arguments name and summaryGraph. See the vignette of this package for details.

Slots

```
id: Object of class "character", the objects unique identifier
name: Object of class "character", the name of the process
type: Object of class "character", the type of process
frameIDs: Object of class "character", the identifiers of the flowSets to which the subprocesses are linkes
summaryGraph: Object of class "qaGraph", a graphical summary of the processe's outcome
frameProcesses: Object of class "list", more detailed information for each flowFrame
```

Methods

```
initialize signature(.Object = "qaProcess"): constructor
```

Author(s)

Florian Hahne

See Also

```
qaGraphList, writeQAReport, qaProcessFrame
```

```
showClass("qaProcess")
```

16 qaProcess.2DStatsPlot

qaProcess.2DStatsPlot $Creates\ a\ QA\ process\ for\ 2D\ summary\ statistic\ comparisons\ across\ multiple\ samples$

Description

This function takes a list of flowSet as input and creates all necessary output for a '2DStatsPlot' type QA process. Objects created by this function can be laid out as HTML using writeQAReport.

Usage

```
qaProcess.2DStatsPlot(flowList,dyes=c("FSC-A","SSC-A"),
  outdir="QAReport",thresh=0.25,func=mean,
  sum.dimensions=NULL,det.dimensions= NULL,pdf=TRUE,
  name="2DStats",...)
```

Arguments

flowList	A list of flowSet
dyes	Flow parameters to be compared across multiple samples
outdir	The directory to which the graphical output is to be saved. If multiple QA processes are to be combined, make sure to use the same directory every time.
thresh	A numeric value between 0 and 1 which is used as a threshold by the outlier detection algorithm
func	Summary statistic function to be applied to each parameter in a flowFrame (mean,median etc)
sum.dimensions	The dimensions of summary plot generated in inches, default is NULL where the dimensions are automatically selected
det.dimensions	the dimensions of each sub image generated by the QA process in inches, default is NULL where the dimensions are automatically selected.
pdf	boolen value determinining if pdf files or jpeg images will be produced by the QA process
name	The name of the process used for the headings in the HTML output.
	Further arguments.

Details

QA processes of type '2DStatsPlot' detects differences in the value of a summary statistic such as the mean, median etc across multiple samples.

For more details on how to layout qaProcess objects to HTML, see writeQAReport and qaReport.

Value

An object of class qaProcess.

Author(s)

Nishant Gopalakrishnan

See Also

 $write QAR eport, \ qaProcess, \ qaProcess. margine vents, \ qaProcess. Density Plot, \ qaProcess. ECDFPlot$

Examples

qaProcess.BoundaryPlot

Creates a QA process for comparison of the percentage of boundary events for the specified parameters across multiple samples

Description

This function takes a list of flowSet as input and creates all necessary output for a 'BoundaryPlot' type QA process. Objects created by this function can be laid out as HTML using writeQAReport.

Usage

```
qaProcess.BoundaryPlot(flowList, dyes=NULL,
  outdir="QAReport", cutoff=3, sum.dimensions=NULL,
  det.dimensions=NULL,pdf=TRUE, name="Boundary",
  side="both",...)
```

Arguments

flowList	A list of flowSet
dyes	Flow parameters to be compared across multiple samples. If set to NULL, all parameters that are duplicated across the list of flowSet are identified and compared.
outdir	The directory to which the graphical output is to be saved. If multiple QA processes are to be combined, make sure to use the same directory every time.
cutoff	A numeric value for the percentage of boundary events that is used by the process to identify samples that failed the QA process
sum.dimensions	The dimensions of summary plot generated in inches, default is NULL where the dimensions are automatically selected
det.dimensions	the dimensions of each sub image generated by the QA process in inches, default is NULL where the dimensions are automatically selected.
pdf	boolen value determinining if pdf files or jpeg images will be produced by the QA process

18 qaProcess.cellnumber

name	The name of the process used for the headings in the HTML output.
side	Object of class "character". The margin on which to evaluate the filter. Either upper for the upper margin or lower for the lower margin or both margins.
	Further arguments.

Details

QA processes of type 'BoundaryPlot' helps identify samples that have a large number of boundary events that must be removed by further gating before detailed analysis of the data.

For more details on how to layout qaProcess objects to HTML, see writeQAReport and qaReport.

Value

An object of class qaProcess.

Author(s)

Nishant Gopalakrishnan

See Also

```
writeQAReport, qaProcess, qaProcess.DensityPlot, qaProcess.ECDFPlot
```

Examples

```
qaProcess.cellnumber Create QA process of type 'cellnumber'
```

Description

This function takes a flowSet as input and creates all necessary output for a 'cellnumber' type QA process. Objects created by this function can be laid out as HTML using writeQAReport.

Usage

qaProcess.cellnumber 19

Arguments

set	A flowSet.
grouping	A character vector defining one of the variables in the phenoData of set used as a grouping variable. If this argument is used and if absolute.value is NULL, outlier detection will be performed within groups rather than across all samples.
outdir	The directory to which the graphical output is to be saved. If multiple QA processes are to be combined, make sure to use the same directory for all of them.
cFactor	The outlier threshold at which the QA criterion is considered to have failed. This is essentially the factor of standard deviations away from the average number of cells per sample, either in both directions if two.sided=TRUE or only towards smaller event numbers if two.sided=FALSE.
absolute.value	An absolute event count below which the QA criterion is considered to be failed. If this argument is not NULL, cFactor and two.sided are ignored.
two.sided	Perform a two-sided outlier detection, i.e., report both unproportionally high and low event numbers.
name	The name of the process used for the headings in the HTML output.
sum.dimensions	The dimensions of the pdf deviced in inches used for the summary plot.
pdf	Logical indicating whether to create vectorized versions of images for this quality process. This should be set to FALSE if disk space is critical, since the pdf versions of the image consume much more space on the hard drive compared to the bitmap version.
• • •	Further arguments.

Details

QA processes of type 'cellnumber' detect aberations in the number of events per sample. These are either determined dynamically as outliers from the typical distribution of event counts for the whole set, or, if absolute.value is not NULL, by an absolute cutoff value. If there is a natural grouping among the samples, this can be specified using the grouping argument. In this case, the outlier detection will be performed within its respective group for a particular sample.

For more details on how to layout qaProcess objects to HTML, see writeQAReport and qaReport.

Value

An object of class qaProcess.

Author(s)

Florian Hahne

See Also

writeQAReport, qaReport, qaProcess, qaProcess.marginevents, qaProcess.timeflow, qaProcess.timeline

```
## Not run:
data(GvHD)
dest <- file.path(tempdir(), "flowQ")
qp <- qaProcess.cellnumber(GvHD, outdir=dest, cFactor=2)</pre>
```

20 qaProcess.DensityPlot

```
qp
## End(Not run)
```

 ${\tt qaProcess.DensityPlot} \quad \textit{Creates a QA process for comparisons of density plots across multiple} \\ samples$

Description

This function takes a list of flowSet as input and creates all necessary outputs for a 'DensityPlot' type QA process. Objects created by this function can be laid out as HTML using writeQAReport.

Usage

Arguments

flowList	A list of flowSet
dyes	Flow parameters to be compared across multiple samples. If set to NULL, all parameters that are duplicated across the list of flowSet are identified and compared.
outdir	The directory to which the graphical output is to be saved. If multiple QA processes are to be combined, make sure to use the same directory every time.
alpha	Outlier mislabeling rate
absolute.value	If provided, the absolute.value is used for the outlier detection and the alpha argument is ignored. defualt is NULL
sum.dimensions	The dimensions of summary plot generated in inches, default is NULL where the dimensions are automatically selected
det.dimensions	the dimensions of each sub image generated by the QA process in inches, default is NULL where the dimensions are automatically selected.
pdf	boolen value determinining if pdf files or jpeg images will be produced by the QA process
name	The name of the process used for the headings in the HTML output.
• • •	Further arguments.

Details

QA processes of type 'DensityPlot' detects differences in the density plots across multiple samples. A summary measure based on the sum of the pairwise KL Distances of density plots is used as a parameter to flag patient panels that are potential outliers.

For more details on how to layout qaProcess objects to HTML, see writeQAReport and qaReport.

Value

An object of class qaProcess.

qaProcess.ECDFPlot 21

Author(s)

Nishant Gopalakrishnan

See Also

writeQAReport, qaProcess, qaProcess.marginevents, qaProcess.ECDFPlot, qaProcess.2DStatsPlot

Examples

qaProcess.ECDFPlot

Creates a QA process for comparisons of ECDF's across multiple samples

Description

This function takes a list of flowSet as input and creates all necessary outputs for an 'ecdfOutliers' type of QA process. Objects created by this function can be laid out as HTML using writeQAReport.

Usage

```
qaProcess.ECDFPlot(flowList, dyes=NULL,outdir="QAReport",
    alpha=0.05, absolute.value=NULL,sum.dimensions=NULL,
    det.dimensions=NULL, pdf=TRUE,
    name="ECDF",...)
```

Arguments

flowList	A list of flowSet
dyes	Flow parameters to be compared across multiple samples. If set to NULL, all parameters that are duplicated across the list of flowSet are identified and compared.
outdir	The directory to which the graphical output is to be saved. If multiple QA processes are to be combined, make sure to use the same directory every time.
alpha	Outlier mislabeling rate
absolute.value	If provided, the absolute value is used for the outlier detection and the alpha argument is ignored. defualt is $NULL$
sum.dimensions	The dimensions of summary plot generated in inches, default is NULL where the dimensions are automatically selected
det.dimensions	the dimensions of each sub image generated by the QA process in inches, default is NULL where the dimensions are automatically selected.

pdf	boolen value determining if pdf files or jpeg images will be produced by the
	QA process
name	The name of the process used for the headings in the HTML output.
	Further arguments.

Details

QA processes of type 'ECDFPLot' detects differences in the empirical cumulative distribution of flow parameters across multiple samples. A summary measure based on the sum of the pairwise KL Distances of ECDF plots across multiple samples is used as a parameter to flag patient panels that are potential outliers.

For more details on how to layout gaProcess objects to HTML, see writeQAReport and gaReport.

Value

An object of class qaProcess.

Author(s)

Nishant Gopalakrishnan

See Also

writeQAReport, qaReport, qaProcess, qaProcess.DensityPlot qaProcess.marginevents, qaProcess.2DStatsPl

Examples

qaProcess.KLDistPlot Creates a QA process for comparisons of KL Distances across multiple samples

Description

This function takes a list of flowSet as input and creates all necessary outputs for an 'KLD-istPlot' type of QA process. Objects created by this function can be laid out as HTML using writeQAReport.

Usage

```
qaProcess.KLDistPlot(flowList,dyes=NULL,outdir="QAReport",
alpha=0.05, absolute.value=NULL,sum.dimensions=NULL,
det.dimensions=NULL, pdf=TRUE,
name="KLDist",...)
```

qaProcess.KLDistPlot 23

Arguments

flowList	A list of flowSet
dyes	Flow parameters to be compared across multiple samples. If set to NULL, all parameters that are duplicated across the list of flowSet are identified and compared.
outdir	The directory to which the graphical output is to be saved. If multiple QA processes are to be combined, make sure to use the same directory every time.
alpha	Outlier mislabeling rate
absolute.value	If provided, the absolute.value is used for the outlier detection and the alpha argument is ignored. defualt is NULL
sum.dimensions	The dimensions of summary plot generated in inches, default is NULL where the dimensions are automatically selected
det.dimensions	the dimensions of each sub image generated by the QA process in inches, default is NULL where the dimensions are automatically selected.
pdf	boolen value determinining if pdf files or jpeg images will be produced by the QA process
name	The name of the process used for the headings in the HTML output.
	Further arguments.

Details

QA processes of type 'KLDistPlot' detects differences in the pairwise KL Distances of the flow parameters across multiple samples.

For more details on how to layout qaProcess objects to HTML, see writeQAReport and qaReport.

Value

An object of class qaProcess.

Author(s)

Nishant Gopalakrishnan

See Also

writeQAReport, qaReport, qaProcess. DensityPlot qaProcess. 2DStatsPlot qaProcess. BoundaryPlot qaProcess.

```
qaProcess.marginevents
```

Create QA process of type 'marginevents'

Description

This function takes a flowSet as input and creates all necessary output for a 'marginevents' type QA process. Objects created by this function can be laid out as HTML using writeQAReport.

Usage

```
qaProcess.marginevents(set, channels=NULL, side ="both", grouping=NULL,
   outdir, cFactor=2, absolute.value=NULL, name="margin events",
   sum.dimensions=NULL, det.dimensions=c(3,1), pdf=TRUE, ...)
```

Arguments

set	A flowSet.
channels	A character vector of channel names for which margin events are to be recorded. Will use all available channels except for the time channel if NULL.
side	Object of class "character". The margin on which to evaluate the filter. Either upper for the upper margin or lower for the lower margin or both for both margins.
grouping	A character vector defining one of the variables in the phenoData of set used as a grouping variable. If this argument is used and if absolute.value is NULL, outlier detection will be performed within groups rather than across all samples.
outdir	The directory to which the graphical output is to be saved. If multiple QA processes are to be combined, make sure to use the same directory for all of them.
cFactor	The outlier threshold at which the QA criterion is considered to have failed. This is essentially the factor of standard deviations away from the average number of cells per sample for the respective channel.
absolute.value	An absolute percentage of margin events below which the QA criterion is considered to be failed. If this argument is not NULL, cFactor is ignored.
name	The name of the process used for the headings in the HTML output.
sum.dimensions	The dimensions of summary plot generated in inches, default is NULL where the dimensions are automatically selected
det.dimensions	the dimensions of each sub image generated by the QA process in inches.
pdf	Logical indicating whether to create vectorized versions of images for this quality process. This should be set to FALSE if disk space is critical, since the pdf versions of the image consume much more space on the hard drive compared to the bitmap version. Setting the argument to FALSE will also speed up the process.
	Further arguments.

qaProcess.timeflow 25

Details

QA processes of type 'marginevents' record the number of events that fall on the margins of the measurement range for each channel. Unproportionally high numbers of such events can indicate problems with the instrument settings or with compensation. These are either determined dynamically as outliers from the typical distribution of margin event percentages for the whole set, or, if absolute.value is not NULL, by an absolute cutoff value. If there is a natural grouping among the samples, this can be specified using the grouping argument. In this case, the outlier detection will be performed within its respective group for a particular sample.

For more details on how to layout qaProcess objects to HTML, see writeQAReport and qaReport.

Value

An object of class gaProcess.

Author(s)

Florian Hahne

See Also

writeQAReport, qaReport, qaProcess, qaProcess.timeline, qaProcess.timeflow, qaProcess.cellnumber

Examples

```
## Not run:
data(GvHD)
dest <- tempdir()
qp <- qaProcess.marginevents(GvHD, channels=c("FSC-H", "SSC-H"),
   outdir=dest)
qp
## End(Not run)</pre>
```

qaProcess.timeflow

Create QA process of type 'timeflow'

Description

This function takes a flowSet as input and creates all necessary output for a 'timeflow' type QA process. Objects created by this function can be laid out as HTML using writeQAReport.

Usage

```
qaProcess.timeflow(set, outdir, cutoff=2, name="time flow",
sum.dimensions=NULL, det.dimensions=c(3,2), pdf=TRUE, ...)
```

26 qaProcess.timeflow

Arguments

set	A flowSet.
outdir	The directory to which the graphical output is to be saved. If multiple QA processes are to be combined, make sure to use the same directory for all of them.
cutoff	The threshold at which the QA criterion is considered to have failed. An absolute value in the timeline deviation score as computed by the timeLinePlot function.
name	The name of the process used for the headings in the HTML output.
sum.dimensions	, det.dimensions
	The dimensions of the pdf device in inches used for the summary and the detailed plots.
pdf	Logical indicating whether to create vectorized versions of images for this quality process. This should be set to FALSE if disk space is critical, since the pdf versions of the image consume much more space on the hard drive compared to the bitmap version.
	Further arguments.

Details

QA processes of type 'timeflow' detect disturbances in the flow of cells over time. These indicate problems with the cell suspension, clogging of the needle or similar issues. If the flow rate is to high, the frequency of measuring cell douplettes increases.

For more details on how to layout qaProcess objects to HTML, see ${\tt writeQAReport}\ and\ qaReport.$

Value

An object of class qaProcess.

Author(s)

Florian Hahne

See Also

writeQAReport, qaProcess, qaProcess.marginevents qaProcess.timeline, qaProcess.cellnumber

```
## Not run:
data(GvHD)
dest <- tempdir()
qp <- qaProcess.timeflow(GvHD, outdir=dest, cutoff=1)
qp
## End(Not run)</pre>
```

qaProcess.timeline 27

rocess.timeline Create QA process of type 'timelin
--

Description

This function takes a flowSet as input and creates all necessary output for a 'timeline' type QA process. Objects created by this function can be laid out as HTML using writeQAReport.

Usage

```
qaProcess.timeline(set, channels=NULL, outdir, cutoff=1,
  name="time line", sum.dimensions=NULL, det.dimensions=c(3,2),
  pdf=TRUE, ...)
```

Arguments

set	A flowSet
channels	A character vector of channel names for which the qaReport is to be produced. Will use all available channels except for the time channel if NULL.
outdir	The directory to which the graphical output is to be saved. If multiple QA processes are to be combined, make sure to use the same directory for all of them.
cutoff	The threshold at which the QA criterion is considered to be failed. An absolute value in the timeline deviation score as computed by the timeLinePlot function.
name	The name of the process used for the headings in the HTML output.
sum.dimensions,	${\tt det.dimensions}$ The dimensions of the pdf device in inches used for the summary and the detailed plots.
pdf	Logical indicating whether to create vectorized versions of images for this quality process. This should be set to FALSE if disk space is critical, since the pdf versions of the image consume much more space on the hard drive compared to the bitmap version.
	Further arguments.

Details

QA processes of type 'timeline' detect unusal patterns in the acquisition of fluorescense and light scatter measurements over time. These are detected dynamically by identifying trends in the signal intensity over time or local changes in the measurement intensities. The underlying hypothesis is that measurement values are acquired randomly, hence there shouldn't be any correlation to time.

For more details on how to layout qaProcess objects to HTML, see writeQAReport and qaReport.

Value

An object of class qaProcess.

28 qaProcessFrame-class

Author(s)

Florian Hahne

See Also

writeQAReport, qaProcess, qaProcess.marginevents qaProcess.timeflow, qaProcess.cellnumber

Examples

```
## Not run:
data(GvHD)
GvHD <- transform(GvHD, "FL1-H"=asinh(`FL1-H`), "FL2-H"=asinh(`FL2-H`))
dest <- tempdir()
qp <- qaProcess.timeline(GvHD, channel="FL1-H", outdir=dest, cutoff=1)
qp
## End(Not run)</pre>
```

```
qaProcessFrame-class Class "qaProcessFrame"
```

Description

Abstraction of subitems within a qQA process

Details

This class bundles graphs and aggregators for a single flowFrame. This allows to create processes with subcomponents, where each item in the frameAggregators and frameGraphs lists corresponds to one subprocess, which can be used, for instance, to create individual plots for each flow channel. For QA processes without subcomponents, these slots would simply not be populated.

Objects from the Class

Objects should be created using the constructor:

qaProcessFrame(frameID, summaryAggregator, summaryGraph, frameAggregators, frameGraphs, details where frameID is the ID of the flowFrame the process is linked to, summaryAggregator is an object inheriting from class qaAggregator which summarizes the outcome, summaryGraph is an object of class qaGraph which is the overview graph of the process for the whole frame, details is a list containing any additional information regarding the QA process, frameAggregators is an object of class aggregatorList and frameGraphs is an object of class qaGraphList. The latter two are the collections of aggregators and graphs for each subprocess. Only frameID and summaryAggregator are mandatory arguments.

Slots

```
id: Object of class "character", a unique ID of the object
frameID: Object of class "character", ID of the flowFrame the process is linked to
summaryAggregator: Object of class "qaAggregator", an aggregator summarizing the output of
the process
```

qaProcessSummary 29

```
summaryGraph: Object of class "qaGraph", a graphical summary of the process frameAggregators: Object of class "aggregatorList" a list of aggregators for the subprocesses frameGraphs: Object of class "qaGraphList" a list of graphical summaries for the subprocesses details: A list for any additional information
```

Methods

```
initialize signature(.Object = "qaProcessFrame"): constructor
```

Author(s)

Florian Hahne

See Also

```
qaGraphList, writeQAReport, qaProcess
```

Examples

```
showClass("qaProcessFrame")
```

qaProcessSummary

Class "qaProcessSummary"

Description

An internal class to represent QA summaries

Objects from the Class

The class is internal and not ment for interactive use.

Slots

```
panels: Object of class "list"
ranges: Object of class "matrix"
summary: Object of class "list"
mapping: Object of class "list"
pnams: Object of class "character"
overallSum: Object of class "matrix"
```

Methods

```
writeLines signature(text = "qaGraphList", con = "file"): HTML output
show signature(object = "qaGraphList"): print object details
```

Author(s)

Florian Hahne

30 qaReport

qaReport	Create HTML report using one or several QA process function(s)

Description

This function combines all graphical output of multiple QA process functions for one flowSet in a single hyperlinked HTML document.

Usage

Arguments

set	A flowSet
qaFunctions	A character vector of the names of QA process functions to be used
outdir	The directory to which the HTML report is to be saved.
argLists	lists of argument lists for each of the QA process functions specified via qaFunctions
grouping	A character scalar indicating a variable in the flowSet's phenoData that is used as a grouping factor in the output.
	Further arguments that are passed on to writeQAReport.

Details

This is a simple convenience function to produce HTML QA reports for a single flowSet given a list of QA process functions. For more fine-grained control use function writeQAReport directly. An entry point to the output of this function can be found at outdir/index.html.

Value

The function is called for it's side effects

Author(s)

Florian Hahne

See Also

```
qaProcess.marginevents, writeQAReport, qaProcess, qaProcess.timeline
```

qData 31

End(Not run)

qData

Data for demonstrating QA processes that compares dye parameters from each patient across multiple aliquots

Description

A list of 8 flowSet. Each flowSet has data from four patients and each element in the list corresponds to an aliquot.

Usage

data(qData)

Format

The object contains a list of 8 objects of class flowSet, each composed of 4 flowFrames. Each flowFrame correponds to a sample from an aliquot.

Details

This qData dataset contains flouresence information for samples from four patients. The sample from each patient was separated into eight aliquots and each aliquot was stained using a different combination of dyes. For each patient, some dyes appear more than once in multiple alquots. These along with the Forward and Side scatter information can be used for Quality control.

rangeAggregator-class Class "rangeAggregator"

Description

Abstraction of a range type of aggregator where possible states are within certain ranges (e.g. percentages)

Objects from the Class

Objects can be created by calls of the form new("rangeAggregator", ...) or using the constructor rangeAggregator(x, min, max, passed), where x, min and max are numeric scalars, with x in the range of [min, max], and passed is a logical scalar.

Slots

```
min: Object of class "numeric", the range minimummax: Object of class "numeric", the range maximumx: Object of class "numeric", the value within the rangepassed: Object of class "logical" indicating whether the process has passed the QA requirements
```

Extends

Class numericAggregator, directly. Class qaAggregator, by class "numericAggregator", distance 2.

Methods

Author(s)

Florian Hahne

See Also

```
qaProcess.marginevents, qaReport, qaProcess, qaProcess.timeline, discreteAggregator, factorAggregator, numericAggregator, stringAggregator, binaryAggregator
```

Examples

```
showClass("rangeAggregator")
```

```
stringAggregator-class
```

Class "stringAggregator"

Description

Abstraction of a string type of aggregator for which possible states are indicated through a textual description

Objects from the Class

Objects can be created by calls of the form new("stringAggregator", ...) or using the constructor stringAggregator(x, passed), where x is a character scalar, and passed is a logical scalar.

Slots

```
x: Object of class "character" which is a textual description of the outcome passed: Object of class "logical" indicating whether the process has passed the QA requirements
```

Extends

```
Class "qaAggregator", directly.
```

Methods

```
show signature(object = "stringAggregator"): print object details
writeLines signature(text = "stringAggregator", con = "file", sep = "missing", useBytes= "missing"
write to HTML file connection
```

txtFormatQAObject 33

Author(s)

Florian Hahne

See Also

qaProcess.marginevents, qaReport, qaProcess, qaProcess.timeline, discreteAggregator, factorAggregator, numericAggregator, binaryAggregator, rangeAggregator

Examples

```
showClass("stringAggregator")
```

txtFormatQAObject

Formats results contained in QA process objects into a data frame with appropriate row and column names

Description

Formats information contained in QA process objects into a data frame which can then be easly written to a tab delimited file. Information contained includes the values used in the outlier detection, information whether the particular parameter passed the outlier test and the summary aggregator information.

Usage

```
txtFormatQAObject(qp)
```

Arguments

qр

A qaProcess object

Value

A data. frame containing the values used in outlier detection, thresholds used, information regarding whether the parameter passed/failed the test as well as the summary aggregator information

Author(s)

Nishant Gopalakrishnan

See Also

```
writeQATextReport
```

34 writeQAReport

validProcess

Validate a QAProcess object

Description

Check for integrety and existence of files specified in a qaProcess object.

Usage

```
validProcess(object)
```

Arguments

object

A qaProcess object

Value

The function is called for it's side effects.

Author(s)

Florian Hahne

See Also

qaProcess.marginevents, qaReport, qaProcess, qaProcess.timeline, writeQAReport

 ${\tt writeQAReport}$

Create HTML report for (lists of) qaProcess objects

Description

This function combines all graphical output of multiple QA processes for one or several flowSets in a single hyperlinked HTML document.

Usage

writeQAReport 35

Arguments

A flowSet or a list of several flowSets. set A list of gaProcess objects or, in the case of multiple flowSets, a list of lists processes of gaProcess objects. See below for further details. globalProcess A list of qaProcess objects generated by multipanel comparison qaProcess functions, when set is a list of flowSets. outdir The directory to which the HTML report is to be saved. Each gaProcess object should have been created in the same directory. A character scalar indicating a variable in the flowSet's phenoData that is used grouping as a grouping factor in the output. A logical indicating whether the output should be on one long page, or split over pagebreaks several pages. A logical indicating whether vectorized versions of the images should be inpdf cluded. Setting this to TRUE implies that the individual QA processes have also

Details

Both the information about graphical output generated as part of a QA process as well as the qunatitative or qualitative results are stored in objects of class qaProcess. The creation of such objects is abstracted in dedicated constructor function and the user should call these functions directly rather than creating qaProcess manually. writeQAReport takes lists of such objects and combines their information in a unified HTML document. A grouping factor can be specified to indicate subgroups of the data. In the case of multiple panels, a list of flowSets can be given to writeQAReport, and the function expects a list of lists of processes, where each process list is specific to one panel.

An entry point to the output of this function can be found at outdir/index.html.

been created with setting pdf=TRUE.

Value

The function is mostly called for it's side effects, that is, the generation of a HTML quality report in outdir. A URL to an entry point for browsing of the report is returned.

Author(s)

Florian Hahne

See Also

qaProcess.marginevents, qaReport, qaProcess, qaProcess.timeline

```
## Not run:
data(GvHD)
GvHD <- transform(GvHD, "FL1-H"=asinh(`FL1-H`), "FL2-H"=asinh(`FL2-H`))
dest <- tempdir()
qp1 <- qaProcess.timeline(GvHD, channel="FL1-H", outdir=dest, cutoff=1)
qp2 <- qaProcess.marginevents(GvHD, channels=c("FL1-H", "FL2-H"),
   outdir=dest, cFactor=4)
url <- writeQAReport(GvHD, processes=list(qp1, qp2), outdir=dest,
    grouping="Patient")
browseURL(url)</pre>
```

36 writeQATextReport

End(Not run)

writeQATextReport

Create tab delimited text files for (lists of) qaProcess objects

Description

This function combines all the results of multiple QA processes for one or several flowSets in a single tab delimited text document.

Usage

Arguments

set A flowSet or a list of several flowSets.

processes A list of qaProcess objects or, in the case of multiple flowSets, a list of lists

of qaProcess objects. See below for further details.

globalProcess A list of qaProcess objects generated by multipanel comparison qaProcess

functions, when set is a list of flowSets.

fileName A character vector for the name of the file.

Details

Qunatitative or qualitative results are stored in objects of class qaProcess. writeQATextReport takes lists of such objects and combines their information in a tab delimited text document. In the case of multiple panels, a list of flowSets can be given to writeQAReport, and the function expects a list of lists of processes, where each process list is specific to one panel.

Value

The function is mostly called for it's side effects, that is, the generation of a tab delimited text file. .

Author(s)

Nishant Gopalakrishnan

See Also

```
qaProcess.marginevents, qaReport, qaProcess, qaProcess.timeline
```

writeQATextReport 37

Index

*Topic IO	validProcess, 34
qaReport, 30	*Topic package
writeQAReport, 34	flowQ-package, 2
writeQATextReport, 36	
*Topic classes	aggregatorList, 28
aggregatorList-class, 3	aggregatorList(aggregatorList-class), 3
binaryAggregator-class,4	aggregatorList-class, 3
discreteAggregator-class, 5	
factorAggregator-class, 7	binaryAggregator, 5, 7, 10, 13, 32, 33
numericAggregator-class, 9	binaryAggregator
outlier-class, 10	(binaryAggregator-class), 4
outlierResult-class, 11	binaryAggregator-class,4
qaAggregator-class, 12	
qaGraph-class, 13	discreteAggregator, 4, 7, 10, 13, 32, 33
qaGraphList-class, 14	discreteAggregator
qaProcess-class, 15	(discreteAggregator-class), 5
qaProcessFrame-class, 28	discreteAggregator-class, 5
qaProcessSummary, 29	
rangeAggregator-class, 31	evaluateProcess, 6
stringAggregator-class, 32	6
*Topic datasets	factorAggregator, 4, 5, 10, 13, 32, 33
qData, 31	factorAggregator
*Topic dynamic	(factorAggregator-class), 7
evaluateProcess, 6	factorAggregator-class, 7
locateDuplicatedParameters, 8	flowCore, 3
normailzeSets, 8	flowFrame, 3, 14, 15, 28
preProcessFlowList, 12	flowQ (flowQ-package), 2
qaProcess.2DStatsPlot, 16	flowQ-package, 2
qaProcess.BoundaryPlot, 17	flowSet, 8, 9, 16–27, 30, 35, 36
qaProcess.cellnumber, 18	flowSets, <i>15</i> , <i>34</i> – <i>36</i>
qaProcess.DensityPlot, 20	initialize aggregator ist-mothed
qaProcess.ECDFPlot, 21	<pre>initialize,aggregatorList-method (aggregatorList-class), 3</pre>
qaProcess.KLDistPlot, 22	initialize,qaGraph-method
qaProcess.marginevents, 24	(qaGraph-class), 13
qaProcess.timeflow, 25	initialize,qaGraphList-method
qaProcess.timeline, 27	(qaGraphList-class), 14
qaReport, 30	initialize,qaProcess-method
txtFormatQAObject, 33	(qaProcess-class), 15
writeQAReport, 34	initialize,qaProcessFrame-method
writeQATextReport, 36	(qaProcessFrame-class), 28
*Topic methods	(401 1 000331 1 anic 101033), 20
outliers-methods, 11	list, <i>3</i> , <i>14</i>
*Topic misc	locateDuplicatedParameters, 8, 9, 12

INDEX 39

names,qaGraph-method(qaGraph-class),13	rangeAggregator
normailzeSets, 8	(rangeAggregator-class), 31
normalizeSets, 8	rangeAggregator-class,31
normalizeSets (normailzeSets), 8	
numericAggregator, 4, 5, 7, 13, 32, 33	show,aggregatorList-method
numericAggregator	(aggregatorList-class), 3
(numericAggregator-class), 9	show,binaryAggregator-method
numericAggregator-class,9	(binaryAggregator-class),4
.11	show,discreteAggregator-method
outlier, 11	(discreteAggregator-class), 5
outlier-class, 10	show, factorAggregator-method
outlierResult-class, 11	(factorAggregator-class), 7
outliers (outliers-methods), 11	show,numericAggregator-method
outliers, flowSet-method	(numericAggregator-class),9
(outliers-methods), 11	show, qaGraph-method (qaGraph-class), 13
outliers-methods, 11	show,qaGraphList-method
proDracoccElowlist 12	(qaGraphList-class), 14
preProcessFlowList, 12	show,qaProcess-method
qaAggregator, 3-5, 7, 9, 28, 32	(gaProcess-class), 15
qaAggregator (qaAggregator-class), 12	show, qaProcessFrame-method
qaAggregator-class, 12	(qaProcessFrame-class), 28
qaGraph, 3, 14, 15, 28	show, qaProcessSummary-method
qaGraph (qaGraph-class), 13	(gaProcessSummary), 29
qaGraph-class, 13	show, rangeAggregator-method
qaGraphList, 14, 15, 28, 29	(rangeAggregator-class), 31
qaGraphList (qaGraphList-class), 14	show,stringAggregator-method
qaGraphList-class, 14	(stringAggregator-class), 32
qaProcess, 3–8, 10, 12–23, 25–30, 32–36	stringAggregator, <i>4</i> , <i>5</i> , <i>7</i> , <i>10</i> , <i>13</i> , <i>32</i>
qaProcess (qaProcess-class), 15	stringAggregator
qaProcess-class, 15	(stringAggregator-class), 32
qaProcess. 2DStatsPlot, 16, 21–23	stringAggregator-class, 32
qaProcess.BoundaryPlot, 17, 23	Sti Ingriggi egator Class, 32
qaProcess.cellnumber, 18, 25, 26, 28	timeLinePlot, 26, 27
qaProcess.DensityPlot, 17, 18, 20, 22, 23	txtFormatQAObject, 33
	extrormate, los jecci, so
qaProcess.ECDFPlot, 17, 18, 21, 21	validProcess, 34
qaProcess.KLDistPlot, 22	vector, 3, 14
qaProcess.marginevents, 4, 5, 7, 10, 13, 15,	700001, 5, 17
17, 19, 21, 22, 24, 26, 28, 30, 32–36	warpSet, 8 , 9
qaProcess.timeflow, 19, 25, 25, 28	writeLines,binaryAggregator,file,missing,missing-metho
qaProcess.timeline, 4, 5, 7, 10, 13, 15, 19,	(binaryAggregator-class), 4
25, 26, 27, 30, 32–36	writeLines,data.frame,file,missing,missing-method
qaProcessFrame, 15	(writeQAReport), 34
qaProcessFrame (qaProcessFrame-class),	writeLines, discreteAggregator, file, missing, missing-met
28	(discreteAggregator-class), 5
qaProcessFrame-class, 28	writeLines, factorAggregator, file, missing, missing-metho
qaProcessSummary, 29	(factorAggregator-class), 7
qaProcessSummary-class	writeLines, numericAggregator, file, missing, missing-meth
(qaProcessSummary), 29	(numericAggregator-class), 9
qaReport, 4, 5, 7, 10, 13, 16–23, 25–28, 30,	writeLines,qaProcessSummary,file,missing,missing-metho
32–36	(qaProcessSummary), 29
qData, 31	writeLines,rangeAggregator,file,missing,missing-method
range/agregator 4 5 7 10 13 33	
rangeAggregator, 4, 5, 7, 10, 13, 33	(rangeAggregator-class), 31

40 INDEX

```
writeLines, stringAggregator, file, missing, missing-method (stringAggregator-class), 32 writeQAReport, 3, 6, 14–30, 34, 34, 35, 36 writeQATextReport, 33, 36, 36
```